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METASYS

protect what you need

INSTALLATION AND OPERATING INSTRUCTIONS FOR OIL-FREE COMPRESSORS

ENGLISH

INSTALLATION AND OPERATING INSTRUCTIONS FOR METASYS COMPRESSORS**Versions META Air and META CAM****Table of Contents**

1.	Introduction.....	19
2.	CE directives conformity.....	19
3.	Warning information and symbols.....	19
4.	Transport and storage.....	21
5.	General Check.....	21
6.	Warranty.....	21
7.	Technical datasheet and dimensions	21
8.	Products drawings and parts list.....	21
9.	Product information.....	22
10.	Where to install the compressor.....	23
11.	Pressure gauges installation.....	23
12.	Electrical connection.....	24
13.	Compressor running.....	24
14.	Maintenance scheduling.....	26
15.	Ordinary maintenance periodical checks.....	27
15.1	Condensate draining.....	27
15.2	Cleaning/replacement suction filter.....	27
15.3	Safety valve check.....	27
15.4	Air line and air dryer filter cartridge replacement.....	27
15.5	Replacement of the non-return valve pad	28
16.	Troubles shooting.....	28
APPENDIX A -	Compressors dimensions and technical datasheet.....	99
APPENDIX B -	Electrical diagrams.....	101
APPENDIX C -	Automatic drain.....	107
APPENDIX D -	Electromagnetic compatibility.....	109
APPENDIX E -	EC Declarations.....	112

For further information contact your dealer or METASYS: technical service at the address customerservice@metasys.com or at the numer +43 (0)512 205420 - 510

We also inform you that it is at METASYS customers disposal a huge technical documentation about compressors visiting our web pages at the address www.metasys.com

1. Introduction

This instruction manual has been written to simplify use and maintenance of the compressor. Strict attention should be paid to the contents of this manual to ensure a correct and economical operation of the compressor and the maximum safety for the operator. We strongly recommend the use of original spare parts which will guarantee the efficiency and service life of the compressor.

METASYS reserves the right to make technical modifications without prior notifications.

2. CE directives conformity

Electric oil free air compressors described in this manual are manufactured in conformity with the following EC directives (see enclosed copy of the document):

2006/42/EC	Machine Directive
2014/35/EU	Low Voltage Device Directive
2014/30/EU	EMC Directive
2014/29/EU	Simple Pressure Vessels Directive
93/42/EEC	Medical Devices Directive

3. Warning information and symbols

In the instructions for assembly and use, as well as the packaging and the product itself, use is made of the following terms or symbols to denote data or information of special importance:

	Information, instructions and warnings for the prevention of damage to health or materialise
	Caution! Dangerous electric voltage!
	Caution! Hot surface
	Caution! Compressor can automatically start
	CE marking of medical device
	Handling mark on package: fragile, handle with care
	Handling mark on package: keep away from water
	Handling mark on package: keep away from sunlight
	Mark on package: recyclable material
	Manufacturer
	Manufacturing date

	Dealer
	Code
	Serial Number
	Read the operation instructions manual before mounting
	Dispose of correctly in accordance with EU Directive 2012/19/EU (WEEE). The unit must be disposed of properly. Within the European Union, the unit must be disposed of in accordance with EU Directive 2012/19/EU (WEEE). If you have any questions about the correct disposal of parts, please contact your dental trade supplier.

4. Transport and storage

The compressor is shipped in cardboard that protects the product from damage during transport.

Caution! Always use the original packaging to secure the compressor in the upright position.

Protect the compressor from humidity and extreme temperatures during transport and storage. A compressor in its original packaging can be stored in a warm, dry and dust free area. Max humidity: 70%. Min. Temperature -10°C, Max. Temperature: +40°C.

5. General Check

Remove compressor from the package and check that there are no evident signs or damage and immediately notify to the carrier.

6. Warranty

The warranty period is 12 months if all mandatory ordinary and extraordinary maintenance checks are carried out for the compressor model.

In special cases, the warranty period can be extended by a further 12 months, up to a maximum of 36.

7. Technical datasheet and dimensions

For technical datasheet, dimensions and weight information please check the Appendix A placed in the end of this manual.

8. Products drawings and parts list

For the detailed drawing and part list of the required model, contact your authorized dealer. Documentation is available in the site www.metasys.com.

9. Product information

9.1 Use for the Intended Purpose

The compressor is intended to be used for generating compressed air required for operating dental units of, for similar dental applications and other applications where oil-free compressed air is required.

Installation in medical care facilities:

In designing and constructing the compressor, allowance has been made for the requirements of medical products where applicable. Accordingly, the unit can be used for installation in medical care facilities. If the unit is installed in medical care facilities, the requirements stipulated in Directive 93/42 EEC IEC 60601-1 as well as the relevant norms must be observed as applied to installation and assembly.

9.2 Use other than that for the Intended Purpose

The compressed air produced by the compressor is unsuitable for operating breathing equipment or similar facilities without additional filters required for the operating area.

- The compressors are designed to be operated in dry, ventilated rooms, ambient temperature +5 to +40 °C.
 - The compressors are designed for non-continuous operation: the ratio of the ON cycle time to the ON + OFF cycle (duty cycle) is 50% for the META Air 24/30 and 70% for the META Air and META CAM from 150 to 1300.
 - Do not expose the compressor to rain. The machine must not be operated in a damp or wet environment. Use is also prohibited in proximity to gases or combustible liquids.
 - Prior to installing the compressor in medical facilities, it must be ensured that the available medium complies with the requirements stipulated for the relevant purpose in each individual case. Observe the particulars given in Appendix A "Compressors dimensions and technical datasheet".
- When installing, classification and conformity rating must be carried out by the manufacturer of the ultimate product.
- When installed in medical care facilities, the electrical equipment needs special precautions regarding EMC and needs to be installed according to EMC information
 - The compressor should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the compressor should be observed to verify normal operations in the configuration in which it will be used
 - Please consider that mobile RF communications equipment can effect the compressor electrical equipment
 - Any other use or use beyond what is specified is deemed to be not for the intended purpose. The manufacturer accepts no liability for damage resulting there from. All risk is borne solely by the operator/user.

9.3 Product Description

The compressor generates an oil-free, dry (versions with air dryers only) and filtered compressed air required for operating units or dental equipment.

ASSEMBLY

10. Where to install the compressor

Room where compressor has to be installed should be large, well ventilated and protected from dust and intense cold; a dusty environment will cause damages and difficulties in operation.

If dust goes inside, it may reach the air filter, causing rapid clogging and part will be deposited over the components thereby preventing heat exchange. It is therefore evident that the cleanliness of the installation location is extremely important for the proper operation of the machine, as this will avoid excessive operation and maintenance costs.

To facilitate maintenance operations and create favorable air circulation compressor must have a good amount of free space surrounding it (Min. 30 cm). The room should be equipped with openings towards the outside placed in proximity of the floor and the ceiling, which will allow the natural circulation of the air. If this is not possible, fans or extractors must be installed.

It is not necessary to provide for special foundations or bases. The machine may be simply placed on a level floor. Compressors fitted on fixed standing tanks should not be secured to the ground. METASYS recommends installing 4 vibration-damping supports.

Climatic operating conditions:

- Temperature: from +5°C to +40°C
- Air relative humidity: max 70%

Note: The reduction of air density in relation with altitude H is directly proportional to the reduction of Outlet Air efficiency of the compressor.

Calculating formula for loss of efficiency (l/min) in relation with altitude :

$$\rho_H = \rho_0 \left(1 - \left(\frac{6,5}{288} \right) \cdot H \right)^{4,255}$$

Consider $\rho_0 = 1,226 \text{ Kg/m}^3$

$H [\text{Km}]$ = Operating altitude of the compressor

$\rho_H [\text{Kg/m}^3]$ = Air density at altitude H

$\rho_0 [\text{Kg/m}^3]$ = Air density at sea level

11. Pressure gauges installation

For safety reasons some models are provided with pressure gauges not installed. To correctly fix them on the pressure regulator and the pressure switch, always use teflon material to avoid any air leak risk.

12. Electrical connection

The electrical supply line must stand the load indicated on the motor rating. Earthing must always be used for installer safety.

Earth connection is necessary. Before connecting the compressor to the control panel, a high sensitivity switch should be installed on the wall. Never connect the earthing wire to the neutral pole. Electric drawing of METASYS compressors are in Appendix B at the end of the manual

If any electric cable or air hose is damaged it must be immediately replaced. Electric cable may not contact hot parts of the compressor, insulation could be damaged!

For THREE-PHASE motors: when installing pay attention to the direction of rotation of the started motor, it must be the same as indicated by the arrow on the pump unit.

13. Compressor running

Connect the electric line by plug and the compressor to the air line connecting it to the outlet placed on the filter (position number 19 - see exploded view).

IEC Power plug


Power plug

General switch

In case of soundproofed cabinet versions compressor feeding plug is the cabinet one, whereas the general switch is placed on the front of the cabinet. This switch feeds both the cabinet ventilation system (permanent on some models and thermostatic controlled on other) and the pressure switch.


Turn the switch of the pressure switch (18) in position "1". Running of the compressor is fully automatic and controlled by the pressure switch, which stops it when pressure in the tank (14) reaches a maximum value, allowing compressor to start again when it descends under a fixed pressure value. Read the receiver pressure value on pressure gauge (21). When compressor receiver is under pressure, operator may regulate operating pressure acting on the pressure regulator situated on the filter (19 – if present): to fix pressure required is necessary to rotate the hand-grip regulator in clockwise direction to increase pressure, in counter clockwise direction to reduce it and read the value on its pressure gauge (20). When delivery air pressure required is reached, push on the handgrip in order to block it.

The META Air and CAM versions above 150 not light are equipped with a membrane dryer consisting of a forced ventilation cooling system, triple filtration with condensate separator and 5 and 0.01 micron filters to ensure maximum air purity, and a fully automatic membrane dryer. Correct use of the compressor and periodic replacement of the filter cartridges allow the membrane dryer to operate completely maintenance-free. For further information on maintenance, please contact

customerservice@metasys.com.

On request, the compressor can be equipped with a tank (2 litres) to collect the condensate collected by the separator filter

The META Air 24 and 30 non-light versions are equipped with an adsorption dryer with integrated condensate separation and 0.01 micron filtration system preceded by a forced ventilation cooling system to ensure maximum air purity. Correct use of the compressor allows the adsorption dryer to operate completely maintenance-free. For further information on maintenance, please contact

customerservice@metasys.com.

On request, the compressor can be equipped with a tank (2 litres) to collect the condensate collected by the separator filter.

In case of emergency , disconnect the compressor from the main supply.

Compressors have hot surfaces, contact may cause burns or fire.

Automatic start: compressors automatically start when minimum pressure is reached and stop when maximum pressure is reached.

If during pressure regulation the handgrip doesn't move, do not force it! It's enough to pull it upside.

MAINTENANCE

14. Maintenance scheduling

Operation	Chapter	Periodicity	Performed by
Switch off the compressor at the end of the use		Daily	User
Check that compressor operation does not exceed max allowed duty cycle	9.2	Daily	User
Release condensate (versions without dryer or automatic drain)	15.1	Weekly	User
Suction filter check/replacement	15.2	Yearly / Check every 500 hours	User/Qualified technician
Safety valve check	15.3	Yearly	User/Qualified technician
Periodical ordinary maintenance.	15.4	Yearly / 1500 hours	User/Qualified technician
Check the cleaning (dust, paper or leaves) and the air exchange of the room for proper ventilation of the compressor		Weekly	User
Check the proper operation of the fans		Weekly	User
Replacement of pre-dryer filter cartridges (dryer versions)	15.4	Yearly	User/Qualified technician
Check tightness of joints, overall device examination	Exploded views on https://www.metasys.com/download-area/	Yearly	Qualified technician
Cleaning/check or replacement of the non return valve pad	15.5	Every 2 years/check every year	User/Qualified technician
Check that the capacitor capacity corresponds to the nominal value ($\pm 5\%$)	Exploded views on https://www.metasys.com/download-area/	Every 2 years	Qualified technician
Periodical extra-ordinary maintenance. Piston rings replacement for META Air 24/30	Exploded views on https://www.metasys.com/download-area/	Every 2 years/1.500 hours	Qualified technician
Periodical extra-ordinary maintenance. Piston rings replacement for META Air and META CAM from 150 to 1300	Exploded views on https://www.metasys.com/download-area/	Every 3.000 hours or 4 years	Qualified technician

15. Ordinary maintenance periodical checks

The M versions are equipped with a membrane dryer consisting of a forced-ventilation cooling system, a double filtration of 5 and 0.01 microns to ensure maximum air purity, and a fully automatic membrane dryer. Correct operation of the compressor and periodic replacement of the filter cartridges allow the membrane dryer to operate completely maintenance-free.

The M versions equipped with an adsorption dryer (Air Purification) require no maintenance.


More information can be found at www.metasys.com

Before performing any maintenance on the compressor make sure that the power supply is switched off. disconnect the plug after having turned in "0" position the switch of the pressure switch (pos. 18 in the exploded view). Make also sure that the air tank (14) is released of any pressure check the pressure gauge (21).

15.1 Condensate draining


Control condensate in the air tank (14). Release it at least once a week, switching off the compressor and reducing line pressure till 1 bar. Place a container under the air receiver or close to the drain system, open the drain valve (27) till the complete condensate draining.

15.2 Cleaning/replacement suction filter


Cleaning the intake air filters (30) placed on the proper filter holder on the top of cylinders, removing the holder cover with its seals (unscrew the butterfly screw). Clean them every month with compressed air or water, replace it if necessary.

In case of soundproofed versions it's important to grant a regular compressor cleaning inside the cabinet.

15.3 Safety valve check

Check the proper safety valve (16) operating at the first compressor running. Pull the ring placed at the top of the safety valve, verifying the correct air exit

Warning! Safety valve must not be used to relief air from the air receiver! Always protect eyes from compressed air using eye glasses

15.4 Air line and air dryer filter cartridge replacement.

Dryer pre-filtration and line filters cartridges (if present) must be replaced according to the maintenance scheduling (Par. 14). Follow the present instructions:

- I. Disconnect the compressor from the main supply
- II. Open the air drain cock or the condensate relief from the air receiver and relief residual pressure contained in the air receiver
- III. Manually unscrew the filter receiver
- IV. Unscrew the cartridge as shown in the picture above and replace it with the new cartridge
- V. Fix again the filter receiver

Carefully place the o-ring on the lip of the filter receiver, replace the o-ring if necessary

Models with air dryer are equipped with a pressure safety valve placed before dryer filters, when it operates filters cartridges replacement is required.

15.5 Replacement of the non return valve pad

Remove the valve closing nut with o'ring, spring and pad using a 22mm wrench. Check that the rubber pad is clean; if there are small metal parts or dust, remove them all and clean the flat work surface or replace the pad and secure it carefully to the spring. When finished, tighten the nut to the valve body.

16. Troubles shooting

PROBLEM	PROBABLE CAUSE	REMEDY
Compressor does not start, or stops and does not start again.	Bad connections. Blown fuse. Overload cut-out switch has tripped.	Check connections, verify standard line tension
	No tension or tension too low	Check connections, verify standard line tension
	Air receiver charged	Open drain valve to expel air. The compressor must restart at the minimum pressure set in the pressure switch.
	Solenoid valve does not empty the delivery pipe	Control the solenoid valve, clean or replace it.
	Electric motor capacitor damaged or not properly fitted	Check the tension at the capacitor, in case replace it
Compressor does not reach set pressure and overheats easily.	Inlet air filter is blocked. NOTE: It is also possible that more air is being required than compressor is capable of delivering.	Replace aspiration filters.

PROBLEM	PROBABLE CAUSE	REMEDY
Compressor does not build pressure or very low performance	Compressor valve damaged	Replace the complete valve block including gaskets.
Compressor does not build pressure or very low performance	Compressor piston ring consumed	Verify compressor performance and replace the piston rings
Air leaking from pressure switch valve when compressor is not running (versions without solenoid valve?).	Faulty non-return valve.	Clean or replace the non return valve
Air leaking from solenoid valve when compressor is not running.	Faulty non-return valve.	Clean or replace the non return valve
Air leaking from solenoid valve when compressor is not running (versions with solenoid valve).	Faulty non-return valve	Drain the compressor from the residual pressure. Clean or replace the non return valve
Air leaking from pressure switch valve when compressor is running (compressors without solenoid valve only)	Faulty pressure switch valve	Clean or replace the pressure switch valve
The compressor does not reach the maximum pressure and/or easily overheats	The suction filter is locked. WARNING: it is also possible that the required air is higher than the compressor capacity	Replace the suction filter
	Damage to intake valves	Replace the complete valve plate
	Wear of piston compression rings	Check the thickness of the rings and replace them
Loss of air from the vent valve of the compressor pressure switch off (for versions without the ethrovalve)	Air leakage from the vent valve of the die-off compressor pressure switch	Drain the compressor from the remaining pressure, clean or replace the non-return valve pad
Air loss from the compressor electrovalve off (for non-lettrovalve versions)	Air leakage from the vent valve of the die-off compressor pressure switch	Drain the compressor from the remaining pressure, clean or replace the non-return valve pad
Loss of air from the vent valve of the working compressor pressure switch	Malfunction of the vent valve of the pressure switch	Clean or replace the vent valve.
Air pressure from regulator does not adjust.	Diaphragm inside regulator body is broken	Replace regulator.
Compressor operating, but no air from outlet.	Inlet air filter blocked. Pressure regulator closed. Drain valve open	Replace oil filler/air filter plug. Turn regulator clockwise to set required pressure. Close drain valve.
Safety valve placed before dryer filters is active	Filters are blocked and over pressure condition in the delivery pipe	Replace the cartridges of the dryer's filters.
Electric motor cuts off during normal operating	Electric motor temperature probe switches off the compressor to protect the moto	Too high temperature: verify general conditions

Thermal switch stops the compressor in standard conditions	Thermal switch is damaged	Replace the thermal switch
	Problem with piston rings	Check rings conditions
	Electric motor damaged	Check if compressor starting is regular, replace the motor.

For further information regarding technical support, please contact customerservice@metasys.com.

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Dimensions – weight – air receiver - noise level

Modello - Model	L (mm)	D (mm)	H (mm)	W (Kg)	Air receiver (L –Gal)	NOISE (db/1m)
META Air 24 Light	430	400	600	29	24-6.34	65
META Air 30 Light	600	410	770	38	40-10.57	65
META Air 24	500	420	600	33	24-6.34	65
META Air 30	710	410	770	42	40-10.57	65
META Air 150 Light	600	410	770	46	40-10.57	66
META Air 250 Light	600	410	810	54	40-10.57	66
META Air 450 Light	1100	600	820	113	90-23.77	69
META Air 150	710	410	770	50	40-10.57	66
META Air 250	710	410	810	58	40-10.57	66
META Air 450	1100	600	820	118	90-23.77	69
META Air 650 Light	1100	600	910	173	200-52.83	72
META Air 650	1450	820	900	183	200-52.83	72
META Air 900 Light	1560	1000	1030	220	270-71.3	75
META Air 900	1560	1000	1030	240	270-71.3	75
META Air 1300 Light	1980	780	1050	330	500-132	78

Silenced compressor: dimensions – weight – air receiver - noise level

Modello - Model	L (mm)	D (mm)	H (mm)	W (Kg-IB)	Air receiver (L –Gal)	NOISE (db/1m)
META Air 24	530	500	800	67-147	24-6.34	50
META Air 150 Light – META Air 150	750	530	915	94/98-207/216	40-10.57	53
META Air 250 Light – META Air 250	750	530	915	102/106-225/233	40-10.57	53
META Air 450 Light	1245	725	1020	220-485	90-23.77	56
META Air 450	1245	725	1020	225-495	90-23.77	56

Technical characteristics

META AIR

Modello Model	 LT/MIN (CFM)	 LT/MIN (CFM)	 BAR PSI	 kW Hp	 VOLT/HZ	Max Current absorption (A)	Electric supply Main switch Curve D (A)	RPM
META Air 24 Light, META Air 30 Light	135 (4,8) 145 (5,1)	182 (6,4) 218 (7,7)	7 101	1,1 1,5	230/50 220/60	6,2 6,8	8 8	1370 1700
META Air 150 Light	170 (6,0) 185 (6,5) 170 (6,0)	240 (8,5) 280 (9,8) 240 (8,5)	7 101	1,5 2	230/50 220/60 400/50 3P	10,2 12 4,3	10 14 6	1360 1760 1470
META Air 250 Light	250 (8,8) 270 (9,5) 250 (8,8)	360 (12,8) 420 (14,8) 360 (12,8)	7 101	2,2 3	230/50 220/60 400/50 3P	14,3 16,0 6	14 16 6,3	1350 1800 1470
META Air 450 Light	500 (17,6) 540 (19,0) 500 (17,6)	720 (25,6) 840 (29,6) 720 (25,6)	7 101	4,4 6	230/50 220/60 400/50 3P	28,6 32 12	32 32 14	1350 1800 1470
META Air 650 Light	750 (26,4)	1080 (38,5)	7 101	6,6 9	400/50 3P	18	32	1470
META Air 900 Light	900 (31,6)	1440 (51,2)	7 101	10 13,6	400/50 3P	20	25	1450
META Air 1300 Light	1350	2160	7 101	15	400/50 3P	30	32	1450
META Air 24, META Air 30	120 (3,8) 129 (4,5)	182 (6,4) 218 (7,7)	7 101	1,1 1,5	230/50 220/60	6,3 6,9	8 8	1370 1700
META Air 150	152 (5,3) 165 (5,8) 152 (5,3)	240 (8,5) 280 (9,8) 240 (8,5)	7 101	1,5 2	230/50 220/60 400/50 3P	10,3 12,1 4,4	10 14 6	1360 1760 1470
META Air 250	225 (7,9) 243 (8,6) 225 (7,9)	360 (12,8) 420 (14,8) 360 (12,8)	7 101	2,2 3	230/50 220/60 400/50 3P	14,4 16,1 6,1	14 16 6,3	1350 1800 1470
META Air 450	450 (15,9) 486 (17,1) 450 (15,9)	720 (25,6) 840 (29,6) 720 (25,6)	7 101	4,4 6	230/50 220/60 400/50 3P	28,8 32,2 12,2	32 32 14	1350 1800 1470
META Air 650	660 (30,2)	1080 (38,5)	7 101	6,6 9	400/50 3P	18,3	20	1470
META Air 900	780 (27,5)	1440 (51,2)	7 101	10 13,6	400/50 3P	20,4	25	1450

META CAM

Modello Model	 LT/MIN (CFM)	 LT/MIN (CFM)	 BAR PSI	 kW Hp	 VOLT/HZ	Max Current absorption (A)	Electric supply Main switch Curve D (A)	RPM
META CAM 150, META CAM with cover	170 (6,0) 185 (6,5) 170 (6,0)	240 (8,5) 280 (9,8) 240 (8,5)	10 145	1,5 2	230/50 220/60 400/50 3P	10,2 12 4,3	10 14 6	1360 1760 1470
META CAM 250, META CAM 250 Light, META CAM 250 with cover	250 (8,8) 270 (9,5) 250 (8,8)	360 (12,8) 420 (14,8) 360 (12,8)	10 145	2,2 3	230/50 220/60 400/50 3P	14,3 16,0 6	14 16 6,3	1350 1800 1470
META CAM 450 META CAM 450 with cover	500 (17,6) 540 (19,0) 500 (17,6)	720 (25,6) 840 (29,6) 720 (25,6)	10 145	4,4 6	230/50 220/60 400/50 3P	28,6 32 12	32 32 14	1350 1800 1470
META CAM 650, META CAM 650 Light	660 (30,2)	1080 (38,5)	10 145	10 13,6	400/50 3P	20,4	25	1450

Appendix B

Single phase pressure switch MDR 1-2, models:


META Air/CAM 150

META Air/CAM 250

Single phase pressure switch
MDR 1-2, models:


META Air 24 Light
 META Air 24
 META Air 30 Light
 META Air 30

Single phase TANDEM, models:


META Air/CAM 450 Light

META Air/CAM 450


Three phase TANDEM, models (400V):


META Air/CAM 450 Light
META Air/CAM 450


META Air/CAM 650 Light, META Air/CAM 650 single phase

META Air/CAM 650 Light, META Air/CAM 65 three phase

META Air/CAM 900 Light, META Air/CAM 900 three phase


META Air/CAM 1300 Light

NOTES

The ground conductor PE was omitted
The schematics are subjected
to change without notice

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OFF (A) – Minimo 10 minuti - Minimum 10 minutes

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To set drain frequency (how often air receiver drainage is required : setting from 5 to 40 mins – for oilless compressors never less than 10 minutes).

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To set the lenght of the air receiver draining(setting from 2 to 20 seconds – for oilless compressors never more than 5 seconds).

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Autodrain setting must be done to remove all the water from the air receiver. Correct setting belongs to operating conditions and ambient humidity.

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Periodically clean the internal filter as follows :

1. Close the drain cock (C) to avoid compressed air exit from the air receiver
2. Repeatedly push the TEST (D) button to discharge residual pressure
3. Unscrew the drain cock tap (E)
4. Remove the filter inside the body of the valve and clean it with compressed air
5. Place the filter in the original position, hardly screw again the tap (E)
6. Slowly open again the drain cock (C) and check that there are no air leaks
7. Push the TEST (D) button to check the correct functioning of the device

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For further information please contact the distributor at the address info@metasys.com

APPENDICE D

Appendix D

IEC 60601-1-2 Table 1 Guidance and manufacturer's declaration - electromagnetic emissions		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the end user should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emission - CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emission IEC 61000-3-2	Class A	
Voltage fluctuation/flicker emission IEC 61000-3-3	Complies	

IEC 60601-1-2 Table 2 Guidance and manufacturer's declaration - electromagnetic immunity			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the end user should assure that it is used in such an environment			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for I/O lines	±2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV line to line ±2 kV lines to earth	±1 kV line to line ±2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 5 cycles ☒5% Ut (>95% dip in Ut) for 5 s	☒5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 5 cycles ☒5% Ut (>95% dip in Ut) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an Uninterruptible Power Supply or Battery
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment

IEC 60601-1-2 Table 4
Guidance and manufacturer's declaration - electromagnetic immunity

The equipment is intended for use in the electromagnetic environment specified below. The customer or the end user should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150KHz to 80MHz	3 Vrms	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance. $d=1,167*\sqrt{P}$ $d=1,167*\sqrt{P}$ 80 MHz to 800 MHz $d=2,333*\sqrt{P}$ 800 MHz to 2,5 GHz
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: at 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b) Over the frequency range 150 KHz to 80 MHz, field strength should be less than 3 V/m.

IEC 60601-1-2 Table 6
Recommended distances between portable and mobile RF communication equipment and the device

The equipment is intended for use in the electromagnetic environment in which radiated RF disturbance are controlled. The customer of the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum power of communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to 80 MHz $d=1,17*\sqrt{P}$ m	80 MHz to 800 MHz $d=1,17*\sqrt{P}$ m	800 MHz to 2,5 GHz $d=2,33*\sqrt{P}$ m

IEC 60601-1-2 Table 6 Recommended distances between portable and mobile RF communication equipment and the device			
0,01	0,117	0,117	0,233
0,1	0,370	0,370	0,740
1	1,17	1,17	2,33
10	3,70	3,70	7,40
100	11,7	11,7	23,3
For transmitters rated at maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer Note 1: at 80 MHz and 800 MHz, the higher frequency range applies Note 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people			

EC DECLARATION OF CONFORMITY WITH EC DIRECTIVES

93-42 EEC (based on Att. II of D.Lgs. 46/97)

WE - MGF srl

Pascoli Street, 15 - 20081 ABBIATEGRASSO (MI) – ITALY

DECLARE UNDER OUR RESPONSABILITY FOR MANUFACTURE AND SUPPLY THE PRODUCTS

DENTAL COMPRESSORS models:

META AIR 150 106000015 (1ph, 50Hz); 106000013 (1ph, 60Hz)	META AIR 450 with cover 106000033 (3ph, 60Hz); 106000029 (1ph, 50Hz)
META AIR 150 Light 106000015 (1ph, 50Hz); 106000017 (1ph, 60Hz)	META AIR 650 106000036 (3ph, 50Hz);
META AIR 150 with cover 106000012 (1ph, 50Hz); 106000014 (1ph, 60Hz)	META AIR 650 Light 106000039 (3ph, 60Hz)
META AIR 150 Light with cover 106000016 (1ph, 50Hz); 106000018 (1ph, 60Hz)	META AIR 900 106000040 (3ph, 50Hz);
META AIR 24 106000001 (1ph, 50Hz); 03020210 (1ph, 60Hz)	META AIR 900 Light -
META AIR 24 Light 106000005 (1ph, 50Hz); 106000006 (1ph, 60Hz)	META AIR 1300 Light 03020043 (3ph, 50Hz); 03020046 (3ph, 60Hz)
META AIR 24 with cover 106000002 (1ph, 50Hz); 106000004 (1ph, 60Hz)	META CAM 150 107000001 (1ph, 50Hz); 107000003 (1ph, 60Hz)
META AIR 250 106000022 (3ph, 50Hz); 106000019 (1ph, 50Hz); 106000021 (1ph, 60Hz)	META CAM 150 with cover 107000002 (3ph, 50Hz)
META AIR 250 Light 106000024 (1ph, 50Hz); 106000026 (1ph, 60Hz)	META CAM 250 107000008 (3ph, 50Hz); 107000004 (1ph, 50Hz); 107000006 (1ph, 60Hz)
META AIR 250 Light with cover 106000025 (1ph, 50Hz);	META CAM 250 Light -
META AIR 250 with cover 106000023 (3ph, 50Hz); 106000020 (1ph, 50Hz);	META CAM 250 with cover 107000005 (1ph, 50Hz); 107000007 (1ph, 60Hz); 107000009 (3ph, 50Hz);
META AIR 30 106000007 (1ph, 50Hz); 106000008 (1ph, 60Hz)	META CAM 450 107000010 (3ph, 50Hz);
META AIR 30 Light 106000009 (1ph, 50Hz); 106000010 (1ph, 60Hz)	META CAM 450 with cover 107000011 (3ph, 50Hz);
META AIR 450 106000030 (3ph, 50Hz); 106000027 (1ph, 50Hz); 106000028 (1ph, 60Hz)	META CAM 650 107000013 (3ph, 50Hz);
META AIR 450 Light 106000034 (3ph, 50Hz)	META CAM 650 Light -
META AIR 450 Light with cover 106000035 (3ph, 50Hz)	

TO WHICH THIS DECLARATION RELATES, ARE IN CONFORMITY
WITH FOLLOWING DIRECTIVES EEC 93/42/EEC

risk class IIa, in accordance with rule 11 of Annex IX of Directive 93/42 / EEC and subsequent amendments (implemented in Italy with Legislative Decree 24 February 1997, n.46, and subsequent amendments), amended by Directive 2007/47 / EC (implemented in Italy with Legislative Decree 25 January 2010, n.37):

- complies with the essential requirements and the provisions of Directive 93/42 / EEC and subsequent amendments and additions. as per the Technical File filed at the Company;
- is manufactured in accordance with the Quality System that meets the requirements set out in Annex II (point 4 excluded) of the aforementioned Legislative Decree, as per Certificate no. 0425-MED-004309-00 issued on 24/02/2021 ICIM - Notified Body n. 0425 - Piazza Don Enrico Mapelli, 75 - 20099 Sesto San Giovanni (MI);
- complies with Directive 2011/65 and Directive 2015/863/ EU of the European Parliament and of the Council of 8 March 2021 (RoHS III), on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- complies with Regulation UE 2017/745 about requirements applicable to legacy device with ref, to Art. 120(3)

Valid from serial number: 07163100

Head of technical file: Eng. Davide Fiani

ISSUED IN ABBIATEGRASSO ON THE 01/12/2023

REV. 02/23

Ing. Davide Fiani
Managing Director

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METASYS Medizintechnik GmbH